

## Anti-WFDC2 Antibody (4M577)

### Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4M577
Purification:	Protein G purified

### Applications

Application:	ELISA
--------------	-------

### Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human WAP four-disulfide core domain Protein 2
Antigen Species:	Human
Uniprot ID:	Q14508
Synonyms:	HE4;EDDM4;WAP5;dj461P17.6;WAP four-disulfide core domain 2
Biology Area:	Tags & Cell Markers

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481