

## Anti-Whitlow linker Antibody-FITC (1M598)

### Product Details

Ig Type:	IgG2b
Reactivity:	Human
Conjugation:	FITC
Clone:	1M598
Purification:	Protein G purified

### Applications

Verified Activity:	Untransfected HEK-293F cells (green line) and transfected scFv-based Anti-CD19 CAR containing a Whitlow/218 linker stable cells (red line) were stained with anti-Whitlow linker FITC-conjugated antibody (2µg/1*10 <sup>6</sup> cells).
Application:	FCM
Recommended	FCM:1:50-1:200.

### Properties

Purity:	>95%
Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Whitlow linker-KLH
Synonyms:	Whitlow/218 Linker (E3U7Q)
Biology Area:	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481