

DS34942424

## Chemical Properties

CAS No. :

Formula: C15H17FN2O

Molecular Weight: 260.31

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DS34942424 is an orally potent analgesic that does not exhibit mu opioid receptor agonist activity.
Targets(IC50)	Others,Opioid Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8416 mL	19.2079 mL	38.4157 mL
5 mM	0.7683 mL	3.8416 mL	7.6831 mL
10 mM	0.3842 mL	1.9208 mL	3.8416 mL
50 mM	0.0768 mL	0.3842 mL	0.7683 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481