


## L-Seryl-amp

## Chemical Properties

CAS No. :	52435-67-3
Formula:	C13H19N6O9P
Molecular Weight:	434.30
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	L-Seryl-amp is a peptide.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3026 mL	11.5128 mL	23.0256 mL
5 mM	0.4605 mL	2.3026 mL	4.6051 mL
10 mM	0.2303 mL	1.1513 mL	2.3026 mL
50 mM	0.0461 mL	0.2303 mL	0.4605 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

- Belrhali H, Yaremchuk A, Tukalo M, Berthet-Colominas C, Rasmussen B, Bösecke P, Diat O, Cusack S. The structural basis for seryl-adenylate and Ap4A synthesis by seryl-tRNA synthetase. *Structure*. 1995 Apr 15;3(4):341-52. PubMed PMID: 7613865.
- Kohda K, Tada M, Hakura A, Kasai H, Kawazoe Y. Formation of 8-hydroxyguanine residues in DNA treated with 4-hydroxyaminoquinoline 1-oxide and its related compounds in the presence of seryl-AMP. *Biochem Biophys Res Commun*. 1987 Dec 31;149(3):1141-8. PubMed PMID: 3122745.
- Chimnarok S, Gravers Jeppesen M, Suzuki T, Nyborg J, Watanabe K. Dual-mode recognition of noncanonical tRNAs(Ser) by seryl-tRNA synthetase in mammalian mitochondria. *EMBO J*. 2005 Oct 5;24(19):3369-79. Epub 2005 Sep 15. PubMed PMID: 16163389; PubMed Central PMCID: PMC1276171.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481