

NGFR/p75NTR Protein, Human, Recombinant

General Information

Synonyms:	Gp80-LNGFR;p75 ICD;TNFRSF16;CD271
Protein Construction:	Lys29-Asn250
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P08138
Molecular Weight:	32~60 kDa (Reducing conditions)

QC Testing

Biological Activity:	ED 50 < 0.4 µg /ml, measured in a neutralization assay using TF-1 cells in the presence of 10.0 ng/ml human b-NGF.
Purity:	> 95% as determined by SDS-PAGE
Endotoxin:	< 0.2 EU/µg of protein as determined by the LAL method.
Formulation:	Lyophilized from a 0.2 µm filtered solution in PBS.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Upon receiving, this product remains stable for up to 6 months at lower than -70°C. Upon reconstitution, the product should be stable for up to 1 week at 4°C or up to 3 months at -20°C. For long term storage it is recommended that a carrier protein (example 0.1% BSA) be added. Avoid repeated freeze-thaw cycles.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

NGF Receptor, also known as Gp80-LNGFR, p75 ICD, CD271 and TNFRSF16, is a type I transmembrane protein belonging to the TNF receptor family. It is expressed by both neuronal and non-neuronal cells. Signaling through NGF Receptor has been shown to regulate gene expression, cell migration and death. A truncated NGF Receptor containing only the extracellular domain has been detected in plasma, amniotic fluid and urine, and acts as a potent NGF antagonist.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481