

ADC Control Human IgG1-McMMAF

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	ADC Control Human IgG1-McMMAF is an antibody-drug conjugate (ADC) comprising the ADC antibody IgG1 Human IgG1 kappa, Isotype Control, and the payload McMMAF.
Targets(IC50)	Antibody-Drug Conjugates (ADCs)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481