


PROTAC SMARCA2/4 degrader-36

Chemical Properties

CAS No. :	3033588-00-7
Formula:	C53H62ClN9O4S
Molecular Weight:	956.635
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PROTACSMARCA2/4 degrader-36 (Compound 29) is a potent dual degrader targeting SMARCA2 and SMARCA4, exhibiting DC50 values of 0.22 nM and 0.85 nM, respectively. Additionally, PROTACSMARCA2/4 degrader-36 demonstrates antiproliferative activity.
Targets(IC50)	Epigenetic Reader Domain,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0453 mL	5.2266 mL	10.4533 mL
5 mM	0.2091 mL	1.0453 mL	2.0907 mL
10 mM	0.1045 mL	0.5227 mL	1.0453 mL
50 mM	0.0209 mL	0.1045 mL	0.2091 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481