


Casein kinase 1δ-IN-23

Chemical Properties

CAS No. :	729579-66-2
Formula:	C ₁₉ H ₁₅ N ₃ O ₃ S
Molecular Weight:	365.406
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Casein kinase1δ-IN-23 (compound 423) is an effective inhibitor of casein kinase1δ. It is applicable in research related to neurodegenerative diseases such as Alzheimer's disease.
Targets(IC50)	Casein Kinase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7367 mL	13.6833 mL	27.3665 mL
5 mM	0.5473 mL	2.7367 mL	5.4733 mL
10 mM	0.2737 mL	1.3683 mL	2.7367 mL
50 mM	0.0547 mL	0.2737 mL	0.5473 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481