

MK-761 TFA

Chemical Properties


CAS No. :

Formula: C₁₅H₂₀F₃N₃O₄

Molecular Weight: 363.332

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	MK-761 TFA is an effective, orally active β 2-adrenergic receptor blocker. It exhibits antihypertensive properties and positive inotropic effects, along with vasodilatory activity.
Targets(IC50)	Adrenergic Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7523 mL	13.7616 mL	27.5232 mL
5 mM	0.5505 mL	2.7523 mL	5.5046 mL
10 mM	0.2752 mL	1.3762 mL	2.7523 mL
50 mM	0.055 mL	0.2752 mL	0.5505 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481