

Cy3-PEG-NH2 (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-PEG-NH2 (MW 1000) is a polyethylene glycol derivative that incorporates Cy3 fluorescent dye, polyethylene glycol (PEG), and an amino group. This compound is commonly utilized in applications such as immunolabeling, nucleic acid labeling, fluorescence microscopy, and flow cytometry. The maximum emission wavelength of Cy3 ranges from approximately 562 to 570 nm. Cy3-PEG-NH2 (MW 1000) is suitable for fluorescence imaging, fluorescent tracing, and fluorescent labeling.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481