


## DDR1/2 inhibitor-3

## Chemical Properties

CAS No. :	2241813-33-0
Formula:	C <sub>31</sub> H <sub>31</sub> F <sub>3</sub> N <sub>6</sub> O
Molecular Weight:	560.613
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	DDR1/2 inhibitor-3 (5n) is an inhibitor of DDR1/2, with IC <sub>50</sub> values of 9.4 nM and 20.4 nM, respectively. It is applicable in anti-inflammatory research.
Targets(IC <sub>50</sub> )	Discoidin Domain Receptor (DDR)

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7838 mL	8.9189 mL	17.8377 mL
5 mM	0.3568 mL	1.7838 mL	3.5675 mL
10 mM	0.1784 mL	0.8919 mL	1.7838 mL
50 mM	0.0357 mL	0.1784 mL	0.3568 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481