


Thalidomide-4-piperidineacetaldehyde

Chemical Properties

CAS No. :	2222115-36-6
Formula:	C ₂₀ H ₂₁ N ₃ O ₅
Molecular Weight:	383.398
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Thalidomide-4-piperidineacetaldehyde is a conjugate of an E3 ubiquitin ligase ligand and linker, utilized in the synthesis of PROTAC H122.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6082 mL	13.0412 mL	26.0824 mL
5 mM	0.5216 mL	2.6082 mL	5.2165 mL
10 mM	0.2608 mL	1.3041 mL	2.6082 mL
50 mM	0.0522 mL	0.2608 mL	0.5216 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481