

AZ14240475

Chemical Properties


CAS No. :

Formula: C₂₃H₁₅ClF₂N₆O₂

Molecular Weight: 480.854

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	AZ14240475 is a potent, selective, brain-penetrant inhibitor of EGFR Ex20Ins mutants (EGFR Ex20Ins mutants) with a pIC ₅₀ of 7.6, playing a significant role in cancer research.
Targets(IC ₅₀)	EGFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0797 mL	10.3983 mL	20.7965 mL
5 mM	0.4159 mL	2.0797 mL	4.1593 mL
10 mM	0.208 mL	1.0398 mL	2.0797 mL
50 mM	0.0416 mL	0.208 mL	0.4159 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481