

## Cy3 MicroRNA Antagomir Negative Control

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cy3 MicroRNA Antagomir Negative Control is a Cy3-labeled MicroRNA Antagomir designed as a negative control. The sequence of this control is based on cel-mir-239b, a sequence commonly used as a control due to its minimal sequence homology with human, mouse, and rat miRNA.
Targets(IC50)	MicroRNA

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481