

## OH-PEG3-C2-Succinamic acid-2-(Methylamino)ethanol

### Chemical Properties

CAS No. :


Formula: C15H30N2O7

Molecular Weight: 350.408

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	OH-PEG3-C2-Succinamic acid-2-(Methylamino)ethanol serves as a PROTAC linker and is utilized in the synthesis of PROTAC Tamoxifen-PEG-Clozapine.
Targets(IC50)	PROTAC Linker

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8538 mL	14.269 mL	28.538 mL
5 mM	0.5708 mL	2.8538 mL	5.7076 mL
10 mM	0.2854 mL	1.4269 mL	2.8538 mL
50 mM	0.0571 mL	0.2854 mL	0.5708 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481