


MTH1 activator-1

Chemical Properties

CAS No. :	2803422-60-6
Formula:	C ₂₉ H ₂₃ F ₃ N ₄ O ₂
Molecular Weight:	516.514
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	MTH1 activator-1 is an MTH1 activator that enhances endogenous MTH1 activity and significantly reduces 8-oxo-dG levels in cellular DNA. It is useful for investigating the upregulation of oxidative damage repair in nucleotide pools and examining biological effects, as well as for studies aiming to delay or prevent tumorigenesis.
Targets(IC50)	DNA/RNA Synthesis

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9361 mL	9.6804 mL	19.3607 mL
5 mM	0.3872 mL	1.9361 mL	3.8721 mL
10 mM	0.1936 mL	0.968 mL	1.9361 mL
50 mM	0.0387 mL	0.1936 mL	0.3872 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481