

GSK-3050002

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	GSK-3050002 (HGS-1035) is a humanized IgG1 antibody that targets and binds to CCL20. It holds potential for inflammation research. For isotype control, refer to Human IgG1 kappa, Isotype Control.
Targets(IC50)	CCR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481