

## SMARCA2/4-ligand-4

## Chemical Properties

CAS No. : 3033592-40-1


Formula: C<sub>22</sub>H<sub>22</sub>ClN<sub>3</sub>O

Molecular Weight: 379.883

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	SMARCA2/4-ligand-4 is a target protein ligand for PROTACs (Ligands for Target Protein for PROTACs). It is utilized in the synthesis of PROTAC SMARCA2/4 degrader-36.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6324 mL	13.1621 mL	26.3241 mL
5 mM	0.5265 mL	2.6324 mL	5.2648 mL
10 mM	0.2632 mL	1.3162 mL	2.6324 mL
50 mM	0.0526 mL	0.2632 mL	0.5265 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481