

ER ligand-6

Chemical Properties

CAS No. : 2504914-79-6


Formula: C₂₉H₃₂FNO₄

Molecular Weight: 477.567

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	ER ligand-6 serves as the target protein ligand for the PROTAC ER Degrader-14.
Targets(IC50)	Estrogen Receptor/ERR,Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0939 mL	10.4697 mL	20.9393 mL
5 mM	0.4188 mL	2.0939 mL	4.1879 mL
10 mM	0.2094 mL	1.047 mL	2.0939 mL
50 mM	0.0419 mL	0.2094 mL	0.4188 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481