

IPS-06061

Chemical Properties

CAS No. :

Formula: C₂₂H₂₆O₃


Molecular Weight: 338.44

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	IPS-06061 is an orally active molecular glue that forms a ternary complex with CRBN and KRASG12D, capable of degrading KRAS G12D, with a DC50 of less than 500 nM.
Targets(IC50)	Molecular Glues,Ras

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9547 mL	14.7737 mL	29.5473 mL
5 mM	0.5909 mL	2.9547 mL	5.9095 mL
10 mM	0.2955 mL	1.4774 mL	2.9547 mL
50 mM	0.0591 mL	0.2955 mL	0.5909 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481