


BLT-4

## Chemical Properties

CAS No. : 251917-79-0  
 Formula: C<sub>18</sub>H<sub>16</sub>N<sub>2</sub>O<sub>2</sub>  
 Molecular Weight: 292.332  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BLT-4 is a specific and reversible inhibitor of scavenger receptor class B type I (SR-BI). It effectively blocks the SR-BI-mediated transfer of lipids between high-density lipoprotein (HDL) and cells.
Targets(IC50)	Liver X Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4208 mL	17.104 mL	34.2079 mL
5 mM	0.6842 mL	3.4208 mL	6.8416 mL
10 mM	0.3421 mL	1.7104 mL	3.4208 mL
50 mM	0.0684 mL	0.3421 mL	0.6842 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481