

## Vensobafusp alfa

## Chemical Properties

CAS No. : 2724922-85-2

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Vensobafusp alfa (KP-104) is a fusion protein composed of an IgG4 monoclonal antibody targeting complement protein C5, combined with the domains 1-5 of complement factor H (FH1-5). It exhibits anti-inflammatory and immunomodulatory properties. The isotype control for Vensobafusp alfa can be referred to as human IgG4 (S228P) kappa.
Targets(IC50)	Complement System

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481