


FK-565

## Chemical Properties

CAS No. : 79335-75-4  
 Formula: C<sub>22</sub>H<sub>38</sub>N<sub>4</sub>O<sub>9</sub>  
 Molecular Weight: 502.559  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	FK-565 is a ligand of nucleotide-binding oligomerization domain-1 (NOD1) that induces a mouse model of arteritis.
Targets(IC50)	NOD-like Receptor (NLR)

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9898 mL	9.9491 mL	19.8981 mL
5 mM	0.398 mL	1.9898 mL	3.9796 mL
10 mM	0.199 mL	0.9949 mL	1.9898 mL
50 mM	0.0398 mL	0.199 mL	0.398 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481