

VRC-01

## Chemical Properties

CAS No. : 1412901-55-3

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	VRC-01 is a human-derived IgG1 antibody that targets the envelope glycoprotein gp120 protein [gp120]. By partially mimicking the interaction between the CD4 receptor and the HIV-1 gp120 envelope glycoprotein, VRC-01 inhibits viral entry. For the isotype control of VRC-01, refer to HumanIgG1kappa, Isotype Control.
Targets(IC50)	HIV Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481