

Cy3-PEG-DBCO (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-PEG-DBCO (MW 5000) is a fluorescent probe containing the Cy3 dye. It features a clickable functional group, DBCO, and PEG, enabling strain-promoted azide-alkyne cycloaddition (SPAAC) reactions with molecules that have azide groups.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481