


## Gardiquimod hydrochloride

## Chemical Properties

CAS No. :	2956183-81-4
Formula:	C <sub>17</sub> H <sub>24</sub> ClN <sub>5</sub> O
Molecular Weight:	349.858
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Gardiquimod (hydrochloride) is an imidazoquinoline class TLR7/8 agonist. It can inhibit HIV-1 infection in macrophages and activated peripheral blood mononuclear cells (PBMCs). At concentrations below 10 $\mu$ M, Gardiquimod (hydrochloride) specifically activates TLR7.
Targets(IC50)	HIV Protease,TLR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8583 mL	14.2914 mL	28.5829 mL
5 mM	0.5717 mL	2.8583 mL	5.7166 mL
10 mM	0.2858 mL	1.4291 mL	2.8583 mL
50 mM	0.0572 mL	0.2858 mL	0.5717 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481