

E3 Ligase Ligand-linker Conjugate 156

Chemical Properties

CAS No. :


Formula: C₁₉H₂₆N₄O₄

Molecular Weight: 374.434

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3 Ligase Ligand-linker Conjugate 156 is a CRBN E3 ligase ligand-linker conjugate. This compound is utilized for synthesizing PROTACXPO1 degrader-1.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6707 mL	13.3536 mL	26.7073 mL
5 mM	0.5341 mL	2.6707 mL	5.3415 mL
10 mM	0.2671 mL	1.3354 mL	2.6707 mL
50 mM	0.0534 mL	0.2671 mL	0.5341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481