

Cy3-PEG-SH (MW 10000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Cy3-PEG-SH (MW 10000) (Cy3-PEG-Thiol) is a near-infrared fluorescent labeling reagent that integrates Cy3 fluorescent dye, polyethylene glycol (PEG), and thiol (SH) polyethylene glycol derivatives, making it suitable for molecular conjugation. It absorbs at a wavelength of 550 nm and offers high solubility and terminal substitution
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481