

CDK2 degrader 4

Chemical Properties


CAS No. : 2924122-01-8

Formula: C₂₃H₂₆ClN₃O₅

Molecular Weight: 459.923

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CDK2 degrader4 (compound 104) is a potent degrader of CDK2, showing promise for cancer research applications.
Targets(IC50)	CDK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1743 mL	10.8715 mL	21.7429 mL
5 mM	0.4349 mL	2.1743 mL	4.3486 mL
10 mM	0.2174 mL	1.0871 mL	2.1743 mL
50 mM	0.0435 mL	0.2174 mL	0.4349 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481