


Acaricidal agent-1

Chemical Properties

CAS No. :	123293-78-7
Formula:	C ₁₀ H ₁₁ NO ₂
Molecular Weight:	177.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Acaricidal agent-1 (Compound 34) exhibits acaricidal activity, inhibiting <i>P. cuniculi</i> with an LC ₅₀ of 0.82 mM (146.6 µg/mL).
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6433 mL	28.2167 mL	56.4334 mL
5 mM	1.1287 mL	5.6433 mL	11.2867 mL
10 mM	0.5643 mL	2.8217 mL	5.6433 mL
50 mM	0.1129 mL	0.5643 mL	1.1287 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481