

GEN-1053

Chemical Properties

CAS No. : 2733575-82-9

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	GEN-1053 (BNT-313) is a humanized IgG1 antibody designed to target CD27. Upon binding, it forms a hexameric structure with CD27. For isotype control, refer to Human IgG1 lambda 1, Isotype Control.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481