


BRN3OMe

## Chemical Properties

CAS No. : 188349-22-6  
 Formula: C<sub>7</sub>H<sub>13</sub>N<sub>3</sub>O<sub>4</sub>  
 Molecular Weight: 203.196  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	BRN3OMe is an effective inhibitor of HIV-1 reverse transcriptase (HIV-1 reverse transcriptase), showing potential for antiviral drug research.
Targets(IC50)	HIV Protease

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9213 mL	24.6063 mL	49.2126 mL
5 mM	0.9843 mL	4.9213 mL	9.8425 mL
10 mM	0.4921 mL	2.4606 mL	4.9213 mL
50 mM	0.0984 mL	0.4921 mL	0.9843 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481