

ER ligand-9


Chemical Properties

CAS No. : 3026011-20-8

Formula: C₃₁H₃₃NO₃

Molecular Weight: 467.599

Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	ER ligand-9 is a conjugate of an estrogen receptor (Estrogen Receptor/ERR) ligand and a linker, utilized in the synthesis of PROTACs ERD-1233.
Targets(IC50)	Estrogen Receptor/ERR, Target Protein Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1386 mL	10.6929 mL	21.3858 mL
5 mM	0.4277 mL	2.1386 mL	4.2772 mL
10 mM	0.2139 mL	1.0693 mL	2.1386 mL
50 mM	0.0428 mL	0.2139 mL	0.4277 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481