

## Cy3-PEG-FA (MW 5000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cy3-PEG-FA (MW 5000) is a fluorescently labeled biomolecule, consisting of Cy3 fluorescent dye, polyethylene glycol (PEG), and folic acid (FA). The Cy3 fluorophore is commonly applied in immunolabeling and nucleic acid labeling. Its maximum emission wavelength ranges from approximately 562 to 570 nm.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481