

MB-66

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	MB-66 (MAPP-66) is a fully human IgG1 antibody targeting HSV and HIV. It is a monoclonal antibody membrane intended for vaginal application. The isotype control for MB-66 can be referred to as HumanIgG1kappa, Isotype Control.
Targets(IC50)	HIV Protease,HSV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481