

## Volrustomig

### Chemical Properties

CAS No. : 2407760-40-9

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Volrustomig is a bi-engineered fragment crystallizable (Fc) domain, a monovalent bispecific IgG1 monoclonal antibody targeting the key immune checkpoint receptors PD-1 and CTLA-4. It boosts T-cell activation and antitumor immunity, making it a promising immunotherapy for various cancers. Molecular weight: 146.77 kDa.
Targets(IC50)	PD-1/PD-L1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481