

M701

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	The compound M701 is a bispecific antibody that targets CD3 on T cells and EpCAM on tumor cells, aiming to stimulate T cell-mediated tumor cell death. It is being developed for the treatment of malignant ascites by redirecting T cells to attack tumors expressing EpCAM. [Molecular weight: 125.84 kDa]
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481