

Emb-06

Chemical Properties

CAS No. : 2919209-65-5

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Emb-06 is a tetravalent CD3 × BCMA bispecific antibody in the form of Fabs-In-Tandem-Ig, which possesses four active antigen-binding sites. It targets BCMA on tumor cells and CD3 on T cells, enhancing immune activation and antitumor efficacy. The molecular weight is 241.36 KD.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481