

Hb0025

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Hb0025 is a bispecific antibody targeting PD-L1 and vascular endothelial growth factor (VEGF), with EC50 values of 0.67 nM and 0.48 nM, respectively. It exerts antitumor effects by inhibiting VEGF activity and facilitating T-cell infiltration into tumors, thereby transforming the immunosuppressive tumor microenvironment into an immune-activated state. Molecular weight: 170.8 KD.
Targets(IC50)	VEGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481