

Regn7075

## Chemical Properties

CAS No. : 2765077-69-6

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Regn7075 is an innovative first-in-class bispecific antibody (bsAb) designed to restore immune sensitivity in traditionally unresponsive tumors by connecting CD28+ T cells with tumor cells expressing EGFR. It enhances immune response and antitumor immunity. The molecular weight is 145.24 KD.
Targets(IC50)	EGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481