

Ngm707

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ngm707 is a humanized monoclonal antibody that acts as a dual antagonist against the key inhibitory immune checkpoint receptors ILT2 and ILT4. It shows potential for immunomodulatory and antitumor effects and can be used alone or in combination to treat advanced solid tumors. Molecular weight: 148.32 kDa.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481