

Regn5678

Chemical Properties

CAS No. : 2657613-60-8

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	REGN5678 is a bispecific antibody that targets both PSMA on prostate cancer cells and CD28 on T cells, enhancing immune activation. This drug is being developed to treat prostate cancer by redirecting T cells to tumors expressing PSMA, thereby promoting tumor cell death. Its molecular weight is 145.86 kDa.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481