

HR97 TFA

Chemical Properties

CAS No. :

Formula: C64H115N27O14S.xC2HF3O2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

HR97 TFA is a cell-penetrating peptide that binds with engineered melanin to produce the eye drop solution HR97-SunitiGel. The peptide-drug conjugate within HR97-SunitiGel allows for sustained ocular drug delivery.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481