


(Rac)-Nanatinostat

Chemical Properties

CAS No. :	914937-68-1
Formula:	C ₂₀ H ₁₉ FN ₆ O ₂
Molecular Weight:	394.402
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(Rac)-Nanatinostat ((Rac)-CHR-3996, example 44) is a potent HDAC inhibitor with an IC ₅₀ of less than 330 nM. It exhibits anticancer properties, effectively inhibiting the growth of HeLa, U937, and HUT cells.
Targets(IC ₅₀)	HDAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5355 mL	12.6775 mL	25.355 mL
5 mM	0.5071 mL	2.5355 mL	5.071 mL
10 mM	0.2535 mL	1.2677 mL	2.5355 mL
50 mM	0.0507 mL	0.2535 mL	0.5071 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481