


## Adustin

## Chemical Properties

CAS No. :	69579-74-4
Formula:	C <sub>11</sub> H <sub>8</sub> O <sub>3</sub>
Molecular Weight:	188.179
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Adustin is an antifungal antibiotic characterized as a polypeptide that inhibits translation. It effectively suppresses translation in a cell-free rabbit reticulocyte lysate system, exhibiting an IC <sub>50</sub> of 0.34 μM.
Targets(IC <sub>50</sub> )	Antibiotic, Antifungal

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3141 mL	26.5703 mL	53.1406 mL
5 mM	1.0628 mL	5.3141 mL	10.6281 mL
10 mM	0.5314 mL	2.657 mL	5.3141 mL
50 mM	0.1063 mL	0.5314 mL	1.0628 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481