


CRBN ligand-84

Chemical Properties

CAS No. :	3065251-76-2
Formula:	C ₁₁ H ₈ BrClFNO ₂
Molecular Weight:	320.542
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	CRBN ligand-84 is a CRBN-type E3 ubiquitin ligase ligand, suitable for use in the development of PROTAC [proteolysis-targeting chimera].
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1197 mL	15.5987 mL	31.1974 mL
5 mM	0.6239 mL	3.1197 mL	6.2395 mL
10 mM	0.312 mL	1.5599 mL	3.1197 mL
50 mM	0.0624 mL	0.312 mL	0.6239 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481