

YH32367

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	YH32367 (ABL105) is a bispecific antibody targeting HER2 and 4-1BB. It induces the secretion of IFN- γ , resulting in the death of tumor cells when co-cultured with hPBMC and HER2-expressing tumor cells. YH32367 demonstrates notable antitumor activity.
Targets(IC50)	EGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481