

## HBpep(RPY)-SP TFA

## Chemical Properties

CAS No. :

Formula: C155H193N45O38S2.xC2HF3O2

Molecular Weight:

Keep away from moisture

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HBpep(RPY)-SP TFA is an optimized phase-separating peptide capable of liquid-liquid phase separation to form conjugated microdroplets over a wide pH range (6.5-9.0). It is highly efficient in the delivery of a wide range of biomolecules and can be used in difficult-to-transfect cells.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481