


## Arisugacin A

## Chemical Properties

CAS No. :	167114-89-8
Formula:	C <sub>28</sub> H <sub>32</sub> O <sub>8</sub>
Molecular Weight:	496.549
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Arisugacin A is a potent, orally active, and selective inhibitor of acetylcholinesterase (AChE).
Targets(IC50)	Cholinesterase (ChE)

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0139 mL	10.0695 mL	20.139 mL
5 mM	0.4028 mL	2.0139 mL	4.0278 mL
10 mM	0.2014 mL	1.0069 mL	2.0139 mL
50 mM	0.0403 mL	0.2014 mL	0.4028 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481