


Asticolorin A

Chemical Properties

CAS No. :	93376-70-6
Formula:	C33H30O7
Molecular Weight:	538.587
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Asticolorin A is a mycotoxic metabolite derived from the fungus <i>Aspergillus versicolor</i> .
Targets(IC50)	Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8567 mL	9.2835 mL	18.567 mL
5 mM	0.3713 mL	1.8567 mL	3.7134 mL
10 mM	0.1857 mL	0.9283 mL	1.8567 mL
50 mM	0.0371 mL	0.1857 mL	0.3713 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481